

4-Nitroacetanilide

Chemicals

- Acetic acid, konz.
- Sulfuric acid, konz.
- Acetanilide
- Nitric acid, konz.
- Sodium chloride
- Ethyl alcohol

Material

- 500 mL- multiple neck flask
- Stirring motor
- Stirrer with locking
- Cooler with cooling tubes
- Thermometer with joint
- Beaker glasses
- Dropping funnel
- Cooling bath
- Heating unit
- Measuring cylinder
- Erlenmeyer flask
- Round bottom flask
- Funnel with folded filter
- Powder funnel
- Suction filter with rubber cuff
- Feeding bottle
- Vacuum connection
- Spatula
- Porcelain stick, porcelain bowl
- Drying oven
- Precision balance

Safety tips

Acetanilide

- H272 H302 H410
- P210 P280 P273 P301+P312
- HAZARD!!



Nitric acid, konz.

- H272 H314
- P210 P280 P301+P330+P331
P304+P340 P309+P310
- HAZARD!!



Ethyl alcohol:

- H225
- P210 P243 P280
- HAZARD!!



Acetic acid, konz.

- H226 H314
- P210 P243 P280 P301+P330+P331
P304+P340 P309+P310
- HAZARD!!



Sulfuric acid, konz.

- H314, H290
- P280, P301+P330+P331,
P305+P351+P338, P309+P310
- HAZARD!!




4-Nitroacetanilide

- H319 H335 H315
- P280 P302+P352 P304+P340
P305+P351+P338 P309+P311
- ATTENTION!!



Reaction Equation:

Nitration/ Electrophilic substitution at aromates:



Acetanilid
M=135,17 g/mol

4-Nitroacetanilid
M=180,16 g/mol

Experimental procedure

- Prepare a suspension of 25 g of acetanilide and 25 mL of glacial acetic acid in a 500 mL multiple neck flask
- Stir vigorously and add drop wise 50 mL of conc. sulfuric acid (in the beginning make sure that the velocity of drops is faster to cause a rise of temperature to 40° C)
- After cooling with a mixture of ice and sodium chloride below 10° C, add drop wise to the clear solution a mixture of 13 mL conc. nitric acid and 8 mL conc. sulfuric acid (carefully produced previously using ice cooling)
- The temperature must stay below 10° C during the drop wise addition
- After removing the cooling mixture, stir for another hour at room temperature
- Then pour the contents of the flask onto 250 g of crushed ice and stir
- Under repeated stirring leave the crystal pulp for 15 minutes and suction-filter sharply
- Suspend the filter cake in a beaker glass two times with 300 mL of water each, then suction-filter and wash acid-free (approx. 6x100mL)
- Suction-filter the raw product as dry as possible
- Then recrystallize with ethyl alcohol
- Cool the filtrate to 20° C and suck off
- Wash the filter cake two times with 15 mL of cold ethyl alcohol each
- Dry the product to mass consistency at 100° C

Waste disposal:

- Dispose of the filtrates in the container for acid wastes

Analysis:

- Calculate the yield of product regarding acetanilide in grammes and percentage of theory

Preparation list

Chemicals:

- Acetic acid, conc. 25 mL
- Sulfuric acid, conc. 58 mL
- Acetanilide 25 g
- Nitric acid conc. 13 mL
- Sodium chloride for the ice mixture
- Ethyl alcohol approx.. 350 mL

Tools:

- 500 mL- multiple neck flask
- Stirring motor
- Stirrer with locking
- Cooler with cooling tubes
- Thermometer with joint
- Beaker glasses
- Dropping funnel
- Cooling bath
- Heating unit
- Measuring cylinder
- Erlenmeyer flask
- Round bottom flask
- Funnel with folded filter
- Powder funnel
- Suction filter with rubber cuff
- Feeding bottle
- Vacuum connection
- Spatula
- Porcelain stick, porcelain bowl
- Drying oven
- Precision balance

