

Quantitation of manganese as permanganate

Chemicals

- Manganese (II)-sulfate-monohydrate
- Sulfuric acid conc.
- Sulfuric acid $c=1$ mol/L
- Phosphoric acid conc.
- Potassium periodate

Material

- • Volumetric flask
- Volumetric pipette
- Graduated pipette
- Burettes
- Beaker glasses
- Heating device
- Thermometer
- UV/VIS-photometer
- Quartz cuvette
- Pipete aidr
- Squeeze bottle
- Plotting paper or personal computer for evaluation
- Precision balance





Safety tips



- Wear safety goggles



- Wear suitable safety gloves

<u>Manganese(II)-sulfate-monohydrate</u> <ul style="list-style-type: none">• H373 H411• P260 P273 P314• CAUTION! 	<u>Phosphoric acid concz.</u> <ul style="list-style-type: none">• H314• P280 P301+P330+P331• P305+P351+P338 P309+P310• HAZARD! 
<u>Sulfuric acid, concz.</u> <ul style="list-style-type: none">• H314, H290• P280, P301+P330+P331, P305+P351+P338, P309+P310• HAZARD! 	<u>Potassium periodate:</u> <ul style="list-style-type: none">• H272 H319 H335 H315• P210 P280 P302+P352 P304+P340• P305+P351+P338 P309+P311• HAZARD! 

Experimental procedure

- Measuring wavelength 525 nm
- Molar extinction coefficient $\epsilon=1800 \text{ L}/(\text{mol}\cdot\text{cm})$
- Linearity of method $E=0.20 \text{ bis } 0.90$
- Cuvette 1cm
- Blank solution de-ionised water
- Molar masses $M(\text{MnSO}_4\cdot\text{H}_2\text{O})=169.02 \text{ g/mol}$, $M(\text{Mn})=54.94 \text{ g/mol}$
- Stock solution $\beta(\text{Mn}^{2+})=10.00 \text{ mg/ml}$

- Planning of a dilution strategy
- Preparation of calibrating solutions as well as measuring them photometrically
- Preparation of a calibration curve and determination of the concentration of the samples

Processing of the samples:

- Sample with contaminated manganese sulfate $w\approx 90\%$ → determination of the total mass
- Dissolve the sample in 100 mL water and approx. 20 mL sulfuric acid 1 mol/L
- Oxidation

Operation instruction for the oxidation of Mn(II)-ions to Mn(VII)-ions:

- Dilute Mn(II)-containing solutions in 100 mL-beaker glasses with 10 mL water each
- Then add 10 mL conc. sulfuric acid, 5 mL conc. phosphoric acid and 0,4 g potassium periodate
- After adding potassium periodate heat the beaker glasses to 80°C and hold for 30 min
- Cool the solutions and fill them in 100 mL-volumetric flasks

Waste disposal:

- Please dispose of all wastes in the containers for halogen containing solutions

Analysis:

- Calculate the mass of manganese in mg of the obtained sample

Preparation list

Chemicals:

- Sulfuric acid, conc. 120 mL
- Sulfuric acid $c= 1\text{ mol/l}$ 40 mL
- Phosphoric acid conc. 60 mL
- Potassium periodate 10 g
- Manganese sulfate-monohydrate
- Sodium chloride

Material:

- Volumetric flask
- Balance dish
- UV/VIS-photometer
- Quartz cuvettes
- Volumetric pipette
- Beaker glasses
- Measuring pipette
- Spatula
- Measuring cylinder
- Burette
- Pipette aid
- Squeeze bottle
- Drying oven
- Precision balance
- Plotting paper or personal computer with analysis software

Preparation of the sample:

- Weigh out 2.900g-3.200g manganese sulfate-monohydrate in a weighing glass and impurify it with approx. 300mg NaCl

Preparation of the Mn(II)-Ions containing stock solution

- Solve 15.3828 g manganese sulfate-monohydrate with water and 100 mL Sulfuric acid $c= 1\text{ mol/L}$, fill it in a 500 mL volumetric flask und fill up with water
- 1 mL= 10.00 mg Mn