

## Quantification of salicylic acid

### Chemicals

- Salicylic acid
- Sodium hydroxide
- Sodium chloride
- Water, dest.

### Materials

- Measuring flask
- Volumetric pipette
- Measuring pipette
- Burettes
- UV/VIS-Photometer
- Quartz cuvettes
- Pipetting aid
- Fine jet washing bottle
- Graph paper or PC (spreadsheet programm)
- Precision balance

### Safety instructions



- wear safety goggles



- wear adequate safety gloves

#### Salicylic acid

- H302, H318
- P280, P305+P351+P338, P309+P310
- HAZARD!!
- 



#### Sodium hydroxide

- H314
- P280 P301+P330+P331 P305+P351+P338 P309+P310
- HAZARD!!



### Experimental procedure

- |                                  |  |
|----------------------------------|--|
| • Measuring wave length          | 298 nm   |
| • Molar extinction coefficient   | $\epsilon=3500 \text{ L}/(\text{mol}\cdot\text{cm})$ |
| • Linearity of method            | $E=0.1 \text{ bis } 0.9$                             |
| • Cuvette                        | 1cm  |
| • Blank solution                 | deionized water                                      |
| • Molar mass salicylic acid      | $M=138.12 \text{ g/mol}$                             |
| • Stock solution                 | solve 250 mg salicylic acid in 30 mL sodium hydroxid |
| • Stability of solutions         | 4 hours  |
| • Quantity of calibration levels | minimum 5  |

- Design of a dilution strategy
- Preparation of the calibration solutions as well as their photometric measurement
- Dilution of the sample and measurement
- Preparation of a calibration curve and determination of the sample concentration

**Waste disposal:**

- Dispose of all wastes in the container for basic solutions

**Analysis:**

- Calculation of the mass of salicylic acid in mg of the given sample

## Preparation list

### Chemicals:

- Salicylic acid                      approx. 0.5 g
- Sodium hydroxide                approx. 10 g
- Sodium chloride                 approx. 100 mg

### Materials:

- Measuring flask
- Weighing scoop
- UV/VIS-Photometer
- Quartz cuvettes
- Volumetric pipettes
- Beaker glasses
- Measuring pipettes
- Spatula
- Measuring cylinder
- Burettes
- Pipetting aid
- Fine jet washing bottle
- Drying oven
- Precision balance
- Graph paper or PC (spreadsheet programm)

Preparation of the sample solution:

- Weigh 100-150mg salicylic acid in 100 mL-beaker glass and solve with 20 mL sodium hydroxid solution w=10%. Transfer with deionized water in 100 mL-measuring flask and fill up